

Product datasheet for **GA216387**

Mouse Taf2 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Taf2 CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor 2 |
| Symbol: | Taf2 |
| Synonyms: | 150kDa; 4732460C16Rik; AI425886; CIF150; TAF2B; TAFII-150; TAFII150 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)