


Product datasheet for **GA216390**

Mouse Ttl1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttl1 CRISPRa kit - CRISPR gene activation of mouse tubulin tyrosine ligase-like 1
Symbol:	Ttl1
Components:	GA216390G1, Ttl1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216390G2, Ttl1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216390G3, Ttl1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q91V51</u>
Synonyms:	6330444E16Rik; AV014541
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)