

Product datasheet for GA216435

Mouse Cep95 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cep95 CRISPRa kit – CRISPR gene activation of mouse centrosomal protein 95
Symbol:	Cep95
Components:	GA216435G1, Cep95 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216435G2, Cep95 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216435G3, Cep95 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BVV7</u>
Synonyms:	4732496G2IRik; AII18346; AI448335; Ccdc45; D330027A14Rik; F630025I20Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

