


Product datasheet for **GA216443**

Mouse Hook3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hook3 CRISPRa kit - CRISPR gene activation of mouse hook microtubule tethering protein 3
Symbol:	Hook3
Components:	GA216443G1, Hook3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216443G2, Hook3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216443G3, Hook3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BUK6</u>
Synonyms:	5830454D03Rik; AI317159; E330005F07Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)