

Product datasheet for **GA216447**

Mouse **Pik3r5** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Pik3r5 CRISPRa kit - CRISPR gene activation of mouse phosphoinositide-3-kinase regulatory subunit 5
Symbol:	Pik3r5
UniProt ID:	Q5SW28
Synonyms:	AV230647; F730038115Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)