

Product datasheet for **GA216449**

Mouse Ddx11 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Ddx11 CRISPRa kit - CRISPR gene activation of mouse DEAD/H (Asp-Glu-Ala-Asp/His) box helicase 11 |
| Symbol: | Ddx11 |
| UniProt ID: | Q6AXC6 |
| Synonyms: | 4732462I11Rik; CHL1; CHLR1; essa15a; KRG2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)