


Product datasheet for **GA216459**

Mouse Tafa1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tafa1 CRISPRa kit - CRISPR gene activation of mouse TAFA chemokine like family member 1
Symbol:	Tafa1
Components:	GA216459G1, Tafa1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216459G2, Tafa1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216459G3, Tafa1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7TPG8</u>
Synonyms:	C630007B19Rik; FAM19A1 TAFA-1; Tafa-1; TAFA1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)