

Product datasheet for GA216472

Mouse H4f16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hist4h4 CRISPRa kit - CRISPR gene activation of mouse histone cluster 4, H4
Symbol:	H4f16
Components:	GA216472G1, H4f16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216472G2, H4f16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216472G3, H4f16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P62806</u>
Synonyms:	B130044J01Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

