

Product datasheet for **GA216482**

Mouse Lrig3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lrig3 CRISPRa kit - CRISPR gene activation of mouse leucine-rich repeats and immunoglobulin-like domains 3
Symbol:	Lrig3
UniProt ID:	Q6P1C6
Synonyms:	9030421L11Rik; 9130004I02Rik; 9430095K15Rik; mKIAA3016
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)