

Product datasheet for **GA216484**

Mouse Klri2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Klri2 CRISPRa kit - CRISPR gene activation of mouse killer cell lectin-like receptor family I member 2
Symbol:	Klri2
Components:	GA216484G1, Klri2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216484G2, Klri2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216484G3, Klri2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5DT36</u>
Synonyms:	A530090P03Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)