

## Product datasheet for **GA216490**

### Mouse P4ha3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P4ha3 CRISPRa kit - CRISPR gene activation of mouse procollagen-proline, 2-oxoglutarate 4-dioxygenase (proline 4-hydroxylase), alpha polypeptide III
Symbol:	P4ha3
Synonyms:	D930031A02Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)