


Product datasheet for **GA216521**

Mouse Tmem88b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem88b CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 88B
Symbol:	Tmem88b
Components:	GA216521G1, Tmem88b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216521G2, Tmem88b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216521G3, Tmem88b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TYP4</u>
Synonyms:	A230069A22Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)