

Product datasheet for **GA216528**

Mouse Dop1a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dop1a CRISPRa kit - CRISPR gene activation of mouse DOP1 leucine zipper like protein A
Symbol:	Dop1a
Components:	GA216528G1, Dop1a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216528G2, Dop1a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216528G3, Dop1a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	B130005I07Rik; BC035275; C130028L17; D9Ert809e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)