

## Product datasheet for **GA216538**

### Mouse Pgap3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pgap3 CRISPRa kit - CRISPR gene activation of mouse post-GPI attachment to proteins 3
Symbol:	Pgap3
UniProt ID:	<a href="#">A2A559</a>
Synonyms:	CAB2; D430035D22Rik; PerId1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)