

Product datasheet for GA216539

Mouse D5Ert579e activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	D5Ert579e CRISPRa kit - CRISPR gene activation of mouse DNA segment, Chr 5, ERATO Doi 579, expressed
Symbol:	D5Ert579e
Components:	GA216539G1, D5Ert579e gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216539G2, D5Ert579e gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216539G3, D5Ert579e gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80U59</u>
Synonyms:	9030221A05Rik; A930018H20Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

