

## Product datasheet for **GA216547**

### Mouse Tafa4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tafa4 CRISPRa kit - CRISPR gene activation of mouse Tafa4 chemokine like family member 4
Symbol:	Tafa4
UniProt ID:	<a href="#">Q7TPG5</a>
Synonyms:	C130034I18Rik; Fam19a4 Tafa-4; Tafa-4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)