


Product datasheet for **GA216561**

Mouse Lingo4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lingo4 CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat and Ig domain containing 4
Symbol:	Lingo4
Components:	GA216561G1, Lingo4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216561G2, Lingo4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216561G3, Lingo4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q149C3</u>
Synonyms:	A530050P17Rik; LERN4; Lrrn6d
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)