

Product datasheet for **GA216581**

Mouse Spata33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Spata33 CRISPRa kit - CRISPR gene activation of mouse spermatogenesis associated 33
Symbol:	Spata33
Components:	GA216581G1, Spata33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216581G2, Spata33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216581G3, Spata33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8C624</u>
Synonyms:	4732415M23Rik; AI847264
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)