

Product datasheet for **GA216590**

Mouse Tnpo3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnpo3 CRISPRa kit - CRISPR gene activation of mouse transportin 3
Symbol:	Tnpo3
Components:	GA216590G1, Tnpo3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216590G2, Tnpo3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216590G3, Tnpo3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P2B1</u>
Synonyms:	5730544L10Rik; C81142; C430013M08Rik; D6Erttd313e; mKIAA4133; Trn-SR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)