

Product datasheet for **GA216604**

Mouse Serac1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Serac1 CRISPRa kit - CRISPR gene activation of mouse serine active site containing 1
Symbol:	Serac1
Components:	GA216604G1, Serac1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216604G2, Serac1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216604G3, Serac1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3U213</u>
Synonyms:	4930511N22Rik; D17Ert141e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)