

Product datasheet for **GA216623**

Mouse Ppfia2 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Ppfia2 CRISPRa kit - CRISPR gene activation of mouse protein tyrosine phosphatase, receptor type, f polypeptide (PTPRF), interacting protein (liprin), alpha 2 |
| Symbol: | Ppfia2 |
| UniProt ID: | <u>Q8BSS9</u> |
| Synonyms: | 5330438O12; B230207K17Rik; E130120L08Rik; mKIAA4112 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)