

Product datasheet for GA216637

Mouse Hsf5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hsf5 CRISPRa kit – CRISPR gene activation of mouse heat shock transcription factor family member 5
Symbol:	Hsf5
Components:	GA216637G1, Hsf5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216637G2, Hsf5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216637G3, Hsf5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5ND04</u>
Synonyms:	Gm739; MLZ-220
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

