

Product datasheet for **GA216644**

Mouse *Togaram1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Togaram1 CRISPRa kit - CRISPR gene activation of mouse TOG array regulator of axonemal microtubules 1
Symbol:	Togaram1
UniProt ID:	Q6A070
Synonyms:	A430041B07Rik; BC062107; C130034K06; mKIAA0423
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)