

Product datasheet for **GA216678**

Mouse Ep300 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ep300 CRISPRa kit - CRISPR gene activation of mouse E1A binding protein p300
Symbol:	Ep300
UniProt ID:	B2RWS6
Synonyms:	A430090G16; A730011L11; KAT3B; p300; p300 HAT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)