


Product datasheet for **GA216751**

Mouse Zfp335 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp335 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 335
Symbol:	Zfp335
Components:	GA216751G1, Zfp335 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216751G2, Zfp335 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216751G3, Zfp335 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2A5K6
Synonyms:	1810045J01Rik; NIF-1; Nif1; Znf335
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)