

Product datasheet for **GA216763**

Mouse Tafa3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tafa3 CRISPRa kit - CRISPR gene activation of mouse Tafa3 chemokine like family member 3
Symbol:	Tafa3
UniProt ID:	Q7TPG6
Synonyms:	7530404M11Rik; Tafa-3; Tafa3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)