

Product datasheet for **GA216780**

Mouse Skint2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Skint2 CRISPRa kit - CRISPR gene activation of mouse selection and upkeep of intraepithelial T cells 2
Symbol:	Skint2
Components:	GA216780G1, Skint2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216780G2, Skint2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216780G3, Skint2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A7XUX6
Synonyms:	A430070M19; B7S3; OTTMUSG00000008540; Skint-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)