

Product datasheet for **GA216786**

Mouse Ttl10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttl10 CRISPRa kit - CRISPR gene activation of mouse tubulin tyrosine ligase-like family, member 10
Symbol:	Ttl10
Components:	GA216786G1, Ttl10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216786G2, Ttl10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216786G3, Ttl10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A4Q9F3</u>
Synonyms:	4833412E22Rik; 4930595O22Rik; AU022799; BC042775; TtlI5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)