

Product datasheet for **GA216791**

Mouse Shisa3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Shisa3 CRISPRa kit - CRISPR gene activation of mouse shisa family member 3
Symbol:	Shisa3
UniProt ID:	<u>Q3UPR0</u>
Synonyms:	D830007B15Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)