

Product datasheet for **GA216800**

Mouse Taok3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Taok3 CRISPRa kit - CRISPR gene activation of mouse TAO kinase 3
Symbol:	Taok3
UniProt ID:	<u>Q8BYC6</u>
Synonyms:	2900006A08Rik; A130052D22; A430105I05Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)