

Product datasheet for **GA216810**

Mouse Gm5111 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gm5111 CRISPRa kit – CRISPR gene activation of mouse predicted gene 5111
Symbol:	Gm5111
Components:	GA216810G1, Gm5111 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216810G2, Gm5111 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216810G3, Gm5111 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	A630038D02; EG330305
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

