

Product datasheet for GA216827

Mouse Zc3h4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zc3h4 CRISPRa kit - CRISPR gene activation of mouse zinc finger CCCH-type containing 4
Symbol:	Zc3h4
Components:	GA216827G1, Zc3h4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA216827G2, Zc3h4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA216827G3, Zc3h4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q6ZPZ3</u>
Synonyms:	Bwq1; Gm768; Kiaa1064-hp; mKIAA1064
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

