

Product datasheet for GA216855

Mouse Trim67 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trim67 CRISPRa kit – CRISPR gene activation of mouse tripartite motif-containing 67
Symbol:	Trim67
Components:	GA216855G1, Trim67 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216855G2, Trim67 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216855G3, Trim67 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q505D9</u>
Synonyms:	D130049O2IRik; TNL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

