

Product datasheet for GA216875

Mouse Pramel39 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	A430089119Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA A430089119 gene
Symbol:	Pramel39
Components:	GA216875G1, Pramel39 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216875G2, Pramel39 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216875G3, Pramel39 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Gm16522
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)