

Product datasheet for **GA216897**

Mouse Tigd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tigd3 CRISPRa kit - CRISPR gene activation of mouse tigger transposable element derived 3
Symbol:	Tigd3
UniProt ID:	Q7TM95
Synonyms:	MGC60957; MGC67569
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)