

Product datasheet for **GA216898**

Mouse **Kcnk18** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Kcnk18 CRISPRa kit - CRISPR gene activation of mouse potassium channel, subfamily K, member 18
Symbol:	Kcnk18
UniProt ID:	Q6VV64
Synonyms:	Gm781; Tresk; Tresk-2; Trik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)