


Product datasheet for **GA216909**

Mouse Kcp activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kcp CRISPRa kit - CRISPR gene activation of mouse kielin/chordin-like protein
Symbol:	Kcp
Components:	GA216909G1, Kcp gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216909G2, Kcp gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216909G3, Kcp gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3U492
Synonyms:	AW060220; CRIM-2; Crim2; Gm793; KCP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)