


Product datasheet for **GA216923**

Mouse Cog3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cog3 CRISPRa kit - CRISPR gene activation of mouse component of oligomeric golgi complex 3
Symbol:	Cog3
Components:	GA216923G1, Cog3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216923G2, Cog3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216923G3, Cog3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	E430004N23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)