

Product datasheet for GA216929

Mouse Supv3l1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Supv3l1 CRISPRa kit – CRISPR gene activation of mouse suppressor of var1, 3-like 1 (S. cerevisiae)
Symbol:	Supv3l1
Components:	GA216929G1, Supv3l1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216929G2, Supv3l1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216929G3, Supv3l1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80YD1</u>
Synonyms:	6330443E10Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

