

Product datasheet for GA216931

Mouse Tmem241 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem241 CRISPRa kit – CRISPR gene activation of mouse transmembrane protein 241
Symbol:	Tmem241
Components:	GA216931G1, Tmem241 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216931G2, Tmem241 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216931G3, Tmem241 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3UME2
Synonyms:	6030446N20Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

