

Product datasheet for **GA216935**

Mouse Pheta2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pheta2 CRISPRa kit - CRISPR gene activation of mouse PH domain containing endocytic trafficking adaptor 2
Symbol:	Pheta2
Components:	GA216935G1, Pheta2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216935G2, Pheta2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216935G3, Pheta2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q14B98
Synonyms:	C920005C14Rik; IPIP27B; Ses2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)