


Product datasheet for **GA216936**

Mouse Tmem220 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem220 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 220
Symbol:	Tmem220
Components:	GA216936G1, Tmem220 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216936G2, Tmem220 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216936G3, Tmem220 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BP07</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)