

Product datasheet for **GA216985**

Mouse Trim46 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trim46 CRISPRa kit - CRISPR gene activation of mouse tripartite motif-containing 46
Symbol:	Trim46
Components:	GA216985G1, Trim46 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA216985G2, Trim46 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA216985G3, Trim46 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7TNM2</u>
Synonyms:	TRIFIC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)