


## Product datasheet for **GA217014**

### Mouse Cnrip1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Cnrip1 CRISPRa kit - CRISPR gene activation of mouse cannabinoid receptor interacting protein 1                                                                                                                                                                                                               |
| Symbol:       | Cnrip1                                                                                                                                                                                                                                                                                                        |
| Components:   | GA217014G1, Cnrip1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217014G2, Cnrip1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA217014G3, Cnrip1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q5M8N0</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | 1500041B16Rik; 3110054C06Rik; 5330437A18Rik; AI854501; C2orf32                                                                                                                                                                                                                                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)