


Product datasheet for **GA217018**

Mouse Tmem102 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem102 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 102
Symbol:	Tmem102
Components:	GA217018G1, Tmem102 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217018G2, Tmem102 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217018G3, Tmem102 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UPR7</u>
Synonyms:	4632403M07Rik; Cbap; Tmem102-ps
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)