

Product datasheet for **GA217025**

Mouse *Kcnh4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnh4</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily H (eag-related), member 4
Symbol:	<i>Kcnh4</i>
Synonyms:	BEC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

