

Product datasheet for **GA217031**

Mouse *Ccdc177* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Ccdc177</i> CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 177
Symbol:	<i>Ccdc177</i>
Components:	GA217031G1, <i>Ccdc177</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217031G2, <i>Ccdc177</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217031G3, <i>Ccdc177</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3UHB8
Synonyms:	Gm1568
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)