

Product datasheet for **GA217054**

Mouse Tmem106c activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem106c CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 106C
Symbol:	Tmem106c
UniProt ID:	<u>Q80VP8</u>
Synonyms:	AI046681; BC046621; D15Ert405e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)