

Product datasheet for **GA217073**

Mouse Tmem232 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem232 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 232
Symbol:	Tmem232
Components:	GA217073G1, Tmem232 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217073G2, Tmem232 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217073G3, Tmem232 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5K6N0</u>
Synonyms:	E130009J12Rik; Gm940; Tes13-L; Tes13-S
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)