

Product datasheet for **GA217087**

Mouse Tmem237 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem237 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 237
Symbol:	Tmem237
Components:	GA217087G1, Tmem237 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA217087G2, Tmem237 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA217087G3, Tmem237 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3V0J1
Synonyms:	AI853305; Als2cr4; Gm972
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)